



TAKE MEASUREMENT - SITTING POSITION

1


Prepare SOZO for sitting position



A



B



Remove the cable

Disconnect the cable from the SOZO hand and foot plates and remove it from the stand.

Remove SOZOtouch

Remove the SOZOtouch hand plate by sliding it forward and lifting up.

Remove SOZOstep

Remove the SOZOstep foot plate by sliding it backward and lifting up.

2

Place SOZO for sitting position

Place the SOZOstep foot plate on the floor. Place the SOZOtouch hand plate on a table across the patient or on a pillow in the patient's lap. Reconnect the cable making sure the screw on the connector is facing up and proceed with the measurement.

⚠ Do not use a metal table or chair

The patient should:

- Sit back in chair with body in balance and shoulders rolled back.
- Distribute body weight evenly on both hips.
- Bend knees at a right angle, plus or minus 10 degrees.

The weight limit for SOZO:

Standing: 375 lbs, Sitting: 750 lbs

You may leave the white cord in this position to perform measurements standing and sitting.

TIP: Be sure that bare skin is not touching bare skin. **Example:** Patient is wearing shorts and their bare legs are touching each other. **Solution:** Put a towel, paper towel or piece of cloth between the legs.

